INSTRUCTIONS FOR USE

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Prescription device

The device is restricted to sale by or on the order of a licensed physician.

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1. INTENDED USE	3
2. INDICATION FOR USE	3
3. DEVICE DESCRIPTION	3
4. RETURN THE DEVICE	3
5. SPECIFICATIONS	4
6. CONTRAINDICATIONS	5
7. WARNINGS	5
8. PRECAUTIONS	6
9. PARTS	7
9.1 Product description	
9.2 Controls	
10. PREPARING THE DEVICE	8
10.1 Cleaning the device	
11. BEFORE USING THE DEVICE	9
11.1 Preparation	
11.2 Clothing	
12. USING THE DEVICE	9
12.1 Checklist	
12.2 Fitting the device	
12.3 Switching ON the device	
12.4 Starting the session	
12.5 Switching off the device	
12.6 Removing the device	
12.7 Troubleshooting	
13. EXPLANATION OF SYMBOLS	13
14. EMC INFORMATION	14
15. CUSTOMER SUPPORT DETAILS	18
16. MANUFACTURER DETAILS	18

GETTING STARTED

- Read the instructions for use (IFU).
- If you have any questions or concerns, please contact your healthcare professional (HCP).

1. INTENDED USE

BLOOM⁴³ is a medical device intended to improve local circulation of the skin through the application of topical heat.

2. INDICATION FOR USE

BLOOM⁴³ is indicated to provide topical heat to the skin of the breast for the temporary increase of local circulation.

3. DEVICE DESCRIPTION

BLOOM⁴³ is a wearable medical device consisting of a preconditioning unit connected to a battery via a cable. The preconditioning unit integrates a user interface and a temperature-controlled heating element. Designed for preoperative use, **BLOOM⁴³** delivers supraphysiological preconditioning to the skin through the application of topical heat, which temporarily improves blood flow circulation.

The battery can be securely attached to clothing or carried in a pocket for convenience. The preconditioning unit is intended to be placed within a bra, positioned against the breast skin (ideally centered over the nipple-areola complex), and worn for approximately 3 hours. It operates in a series of 30-minute active heating cycles within the range of 40°C (104°F) to 44°C (111.2°F), alternating with non-active heating periods.

4. RETURN THE DEVICE

Do not dispose of the device. After use, the device should be returned, please refer to the instruction provided in the case.

5. SPECIFICATIONS

- Product Name: **BLOOM⁴³**
- Model Number: **B4301**
- Dimensions device: 138 mm x 175 mm x 50 mm
- Weight device: 435g
- Output power: 15W maximum (heating)
- Input power: 10W (5V 2A)
- Approximate skin temperature during active controlled heating:
from 40°C (104°F) to 44°C (111.2°F)
- One full session consists of six cycles of:
 - 30 minutes active controlled heating
 - 30 minutes of non-heating(see Figure 5.1)


Figure 5.1

- Environmental conditions:
 - **Operating condition:** Temperature: +15°C (59°F) to +35°C (95°F)
Relative Humidity: 30% to 75%
Atmospheric Pressure: 700 hPa to 1060 hPa
 - **Storage condition:** Temperature: -20°C (-4°F) to +35°C (95°F)
Relative Humidity: 25% to 70%
Atmospheric Pressure: 500 hPa to 1060 hPa
 - **Transport condition:** Temperature: -20°C (-4°F) to +35°C (95°F)
Relative Humidity: 25% to 75%
Atmospheric Pressure: 500 hPa to 1060 hPa

Allow 30 minutes for the device to reach operating temperature after removal from minimum storage temperature and placement at room temperature (20 °C).

6. CONTRAINDICATIONS

The device should **NOT** be used in the following circumstances:

- Patients with an open wound or compromised skin (e.g., sunburn, rash etc.) unless advised otherwise by a healthcare professional.
- Patients with infected or inflamed skin;
- Patients who are unconscious or incapacitated;
- Patients with an active implantable device (e.g. cardiac pacemaker, defibrillator etc.)

7. WARNINGS

- Do not use the device if it has not been prescribed by your healthcare professional;
- Do not use the device for more than one treatment session;
- Do not use the device post-operatively following breast surgery;
- Do not apply the device directly to broken skin, unless advised otherwise by your healthcare professional;
- Do not use the device, if the package or any part of the device is damaged. Contact customer support at support@plexaa.com;
- Do not charge the device;
- Do not use the device when showering or bathing;
- Do not use the device if it has been immersed in any type of liquid;
- Do not use the device without proper cleaning if it has been used by another individual. Only use the 70% Isopropyl alcohol wipes provided by the manufacturer, or equivalent alcohol wipes to clean the device;
- Do not apply the device to any other part of the body other than the breast;
- Do not wear or use the device when undergoing radiological investigations such as MRI and CT scans;


- Do not use the device if it has been dropped and it is visibly damaged;
- Do not use the device over areas with piercing or jewelry;
- The Preconditioning unit, Cable and Battery are not user-serviceable. Attempting repairs may affect the device's safety and performance.
 - Do not modify, alter or tamper with the device.
 - Do not use the device if it has been modified, altered or tampered with.
- Do not charge or connect other device to the USB-C charging port.
- Any communications equipment of portable RF module and electro-magnetic, such as antenna or mobile, shall be used not closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. To prevent the downgrade performance happened.
- Do not use the device with oils, liniments, salves, ointments or lotions as this may cause skin burns.

8. PRECAUTIONS


- Use of the device is only for its indicated uses as described in this manual; After reading the manual, if you are still unsure about anything, please contact your healthcare professional.
- Only use the device the night before your breast surgery;
- If you experience pain or if your skin develops a skin reaction, blisters or burns, immediately stop using the device and seek medical help.
- The device contains a cable that may cause strangulation or tripping if used improperly. Be careful of the cable when using the device;
- Keep out of reach of children and animals;
- The device should be handled with care as excessive force may result in damage;
- Keep the device away from fire or extreme heat (including radiators);
- Keep the device free of debris , dust and direct sunlight;
- Keep the device away from moisture from a nebulizer or steam kettle;
- Degraded sensors will result of Error signal, please refer to the section 12.7 troubleshooting.

9. PARTS

9.1 Product description

9.2 Controls

Note:

- The illustrations are presented from the perspective of the user looking directly down at it.
- If any LEDs turn yellow showing an error, please refer to the section 11.7 trouble shooting.

10. PREPARING THE DEVICE

10.1 Cleaning the device

The device should be thoroughly cleaned before each use.

Clean the device surfaces thoroughly using the alcohol wipe provided (70% Isopropyl alcohol wipes), paying particular attention to the **area in contact with the skin** (See figure 10.1).


Figure 10.1

Note: Do not immerse the device in water or any other type of liquid.

11. BEFORE USING THE DEVICE

- The device should be used for breast Cup size from B to E.
- The device is intended to be used by breast surgery patients on the night before their scheduled surgery. The users are typically female or male individuals aged 18 years or older.
- The device is provided fully charged.
- Do **NOT** turn "ON" and press "PLAY" before the device is fully in place.
- The device will only allow one full treatment per user.
- The total duration of one full treatment is approximately **3 hours**.

11.1 Preparation

Ensure the breast is clean and dry before applying the device. Do **NOT** apply any moisturizers, including creams and emollients, before using the device.

11.2 Clothing

We recommend that you wear a well-fitting, non-underwired bra, or as per the advice of your healthcare professional.

12. USING THE DEVICE

12.1 Checklist

- BLOOM⁴³ device.
- Well-fitting, non-underwired bra.
- Breast skin is clean.

12.2 Fitting the device

1. Place the preconditioning unit into the bra (from above), aiming to position the nipple in the center of the cup (Figure 12.2.1).


Figure 12.2.1

2. Reposition the bra, if required, to ensure the device is well-supported and in a comfortable position. The correct position of the preconditioning unit is with the controller and cable in the 12 o'clock position (Figure 12.2.2).
3. You can place the battery in your pocket or attach to your clothing using the hook (Figure 12.2.3).

Make sure that the cable is not trapped between the preconditioning unit and breast skin.


Figure 12.2.2



Figure 12.2.3

12.3 Switching ON the device

Important: Before starting the session, please ensure that the device is fitted according to the instructions (see section 12.2).


1. Turn the device **ON** by pressing and holding the power button until all the cycle LEDs are ON. Two short vibrations will indicate that the device has been switched on.

The Status light will then turn green and the cycle lights pulse white. This indicates the device is ready to use.

12.4 Starting the session

1. Start the treatment session by **pressing the PLAY button** on the device.


cycle 1 	A. One short vibration and the first pulsing light will indicate that the first cycle has started, during which the first active heating cycle will occur.
cycle 2 	B. One short vibration and the second pulsing light will indicate that the second cycle has started, during which there will be no heating.
cycle 3 	C. One short vibration and the third pulsing light will indicate that the third cycle has started, during which another active heating cycle will occur.
cycle 4 	D. One short vibration and the fourth pulsing light will indicate that the fourth cycle has started, during which there will be no heating.
cycle 5 	E. One short vibration and the fifth pulsing light will indicate that the fifth cycle has started, during which another active heating cycle will occur.
cycle 6 	F. One short vibration and the sixth (final) pulsing light will indicate that the sixth cycle has started, during which there will be no heating.
finish 	G. Three short vibrations will indicate the sixth cycle has finished ; this indicates that the session is now complete. All of the lights will be ON and NOT pulsing.

Note: If the device is switched off during the session, please refer to the section 12.7 troubleshooting.

12.5 Switching off the device

When the session has been completed, turn the device **OFF** by **pressing and holding down the power button** until the LED lights are off. There will be two short vibrations.

If you do not turn OFF the device, please note that it will automatically turn itself OFF one hour after the session has been completed.

12.6 Removing the device

Carefully remove the device from the bra. Place the device back in its case.

12.7 Troubleshooting

- 1. If the device is switched off during the session and then restarted, it may perform a holding period before starting the heat cycle.
- 2. If you encounter the following Error signals during the session, please remove the device and contact PLEXā customer support.

Error signals:

- Error or Status LED **lights are yellow with five long vibrations** (Figure 12.7.1).
- Error LED **light alone is yellow** (Figure 12.7.2).


Figure 12.7.1


Figure 12.7.2

13. EXPLANATION OF SYMBOLS

	Power Button: ON/OFF
	START
	Manufacturer
	Date of manufacture
	Medical device
	Serial number
	Model number
	Refer to instruction manual/booklet
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive (WEEE)
	Type BF applied parts
Rx Only	Prescription Device only
	Direct Current
IP22	Ingress protection (protection against solids up to 12.5 mm and dripping water when tilted up to 15°)
	Temperature Limit
	Relative Humidity limitation
	Non-ionizing electromagnetic radiation

14. EMC INFORMATION

This equipment has been tested and approved with the requirements of applicable EMC standards IEC/EN 60601-1-2. Medical electrical equipment is limited to provide precaution and protection against the harmful interference in a typical medical use. If this equipment needs to be applied to other electrical equipment, normal operation must be checked before use.

Guidance and Manufacturer’s Declaration – Electrical Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishment, including domestic establishments and those directly connect to the public low-voltage power supply network.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC 61000-3-3	Complies	

- WARNING:**
- 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
 - 2. Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
 - 3. Portable RF communications equipment (including peripheral such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any parts of the **BLOOM⁴³**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer’s Declaration – Electrical
Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	Contact: ±8 kV contact Air: ±2, 4, 8kV air	Contact: Not applicable Air: ±2, 4, 8 kVair	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable Not applicable	Not applicable as internally powered device
Surge IEC 61000-4-5	±0.5, 1 kV (line to line) ±0.5, 1, 2 kV (line to earth)	Not applicable	Not applicable as internally powered device
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	Voltage Dips: 0% U_T ; 0.5 cycle 0% U_T ; 1 cycle 70% U_T ; 25/30 cycles Voltage Interruptions: 0% U_T ; 250/300 cycles	Not applicable	Not applicable as internally powered device
NOTE: U_T is the a.c. mains voltage prior to the application of the test level.			

Conducted RF IEC 61000-4-6	3 Vrms: 0.15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz	3 Vrms: 0.15 MHz – 80 MHz 6 Vrms: Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM Fields IEC61000-4-3	10 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	10V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz No loss of safety or essential performance 3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz No change in operation	Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

15. CUSTOMER SUPPORT DETAILS

If you have any questions about **BLOOM⁴³** during use, please contact customer support at support@plexaa.com.

16. MANUFACTURER DETAILS

PLEX $\overline{aa}$ Ltd registered address:

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Domicilio social de PLEXAA Ltd: